

NP Medical
now offers tiered pricing
with the Hospital Care line of
pressure activated valves.

Reveal value in every valve.



HOSPITAL CARE

CRITICAL CARE

VALUE STREAM

Material

Molding

Assembly

Verification Level

Verification Level

Verification Level

You asked, NP Medical delivers. Introducing the value-priced Hospital Care line of pressure activated valves to meet your set of needs. These quality driven products share the same design and materials as our acclaimed Critical Care valves. We can now offer them at a lower price by eliminating manufacturing verification costs and passing those savings on to you.



We make quality patient care possible

CRITICAL CARE

HOSPITAL CARE

Inlet	Medical Grade Polycarbonate, Gamma Stable			
Outlet	Medical Grade Polycarbonate, Gamma Stable			
Diaphragm	Medical Grade Silicone	Blue	Medical Grade Silicone	Green
CHARACTERISTIC	SPECIFICATION	VERIFICATION LEVEL	SPECIFICATION	VERIFICATION LEVEL
Cracking Pressure	100 - 350 mbar air	+	70 - 410 mbar air	●
High Back Pressure	up to 22.7 bar for 5 min	+	up to 21 bar for 20 s	●
Counterflow Pressure	Validated to ISO 8536-12, Section 6.5	+	Validated to ISO 8536-12, Section 6.5	●
Visual Requirements	cracks or holes	+	cracks or holes	●
	incomplete assembly	+	incomplete assembly	●
	loose particulate > 0.20 mm ²	+	loose particulate > 0.20 mm ²	●
	burns, gouges, sinks, short shots	+	burns, gouges, sinks, short shots	●
	embedded particulate > 0.20 mm ²	+	embedded particulate > 0.20 mm ²	●
	discolorations, scratches, scuffs	+	discolorations, scratches, scuffs	●

NOTE: Available in Luer Lock and Tube Fitments.

+ for the Critical Care Line denotes higher verification levels.



We make quality patient care possible

101 Union Street
 Clinton, MA 01510
 U.S.A.

www.npmedical.com