


# Luer-Activated Valves


## Needle-free access devices

### Features


- Bi-directional flow upon actuation
- Lipid resistance
- Low priming volume
- Clear housing
- 100% in-line tested
- Back pressure rated to 45 psi (3.1 bar)
- Latex and DEHP free



FLL/MLL



FLL/Weld Fitment



FLL/T-Style Tubing Port

Inlet/Outlet	Part Number	Housing Material	Avg. Opening Pressure		Priming Volume
			(psi)	(mbar)	(ml)
FLL/MLL	5N02-D300-SSM41G	PC	-	-	0.14
FLL/WF	5N02-D200-SSM41G	PC	-	-	0.19
FLL/T-style TP; ID 0.149"-0.163" (4mm)	5N02-C500-EEM45G	COPE	-	-	0.30

InventIV. CollaboratIV. TransformatIV.

