

The image shows a blue plastic secondary set bag hanger, which is a long, thin, curved device with a hook at one end and a handle at the other. It is designed to hold a secondary set bag for medical use. The hanger is shown in a close-up view, highlighting its ergonomic design and the hook mechanism. The background is a solid blue color.

Medi-Hook Plastic Secondary Set Bag Hanger

The logo for NP Medical, featuring a stylized blue and black icon resembling a medical symbol or a letter 'N' with a curved element.

NPMedical

We Make Quality Patient Care Possible™

The Medi-Hook has been designed and manufactured to perform a crucial function in IV medication delivery. As a helpful accessory for secondary IV medication sets, the Medi-Hook properly positions the primary solution bag to the appropriate height. It features a compact and ergonomic design to simplify automated kit assembly as well as easy clinical manipulation and set up. Manufactured from lightweight, polyethylene material, the Medi-Hook contributes minimal weight to the packaged kits.



Easy to setup and hang up with large handling loop



Robust design, can withstand 20 lbs. pull force



Compact design allows for automated kit assembly



Folds for easy packaging



Occupies minimal weight and space in assembled kits

SPECIFICATIONS

- **Material:** Polyethylene
- **Pull Test:** Withstands 20 lbs. minimum pull test

DIMENSIONS

- **Folded Length:** 5.22" (132.5mm)
- **Width:** 1.28" (32.5mm)
- **Extended Length:** 7.94" (201.7mm)

MANUFACTURING

- ISO 13485 Certified
- ISO Class 8 Cleanroom Manufacturing

COMMERCIAL AVAILABILITY

- Bulk, non-sterile
- 2,000 unit per case

PART NUMBER

5N02-A300H