


The image features a blue plastic hanger for a secondary set bag, with a hook and a handle. The hanger is shown in two views: a top-down view and a side view. The top-down view shows the hook and the handle with a circular opening. The side view shows the hanger's profile and the hook's shape. The background is a solid blue color.

# Medi-Hook Plastic Secondary Set Bag Hanger



We Make Quality Patient Care Possible™

The Medi-Hook has been designed and manufactured to perform a crucial function in IV medication delivery. As a helpful accessory for secondary IV medication sets, the Medi-Hook properly positions the primary solution bag to the appropriate height. It features a compact and ergonomic design to simplify automated kit assembly as well as easy clinical manipulation and set up. Manufactured from lightweight, polyethylene material, the Medi-Hook contributes minimal weight to the packaged kits.



Easy to setup and hang up with large handling loop



Robust design, can withstand 20 lbs. pull force



Compact design allows for automated kit assembly



Folds for easy packaging



Occupies minimal weight and space in assembled kits

### SPECIFICATIONS

- **Material:** Polyethylene
- **Pull Test:** Withstands 20 lbs. minimum pull test

### DIMENSIONS

- **Folded Length:** 5.22" (132.5mm)
- **Width:** 1.28" (32.5mm)
- **Extended Length:** 7.94" (201.7mm)

### MANUFACTURING

- ISO 13485 Certified
- ISO Class 8 Cleanroom Manufacturing

### COMMERCIAL AVAILABILITY

- Bulk, non-sterile
- 2,000 unit per case

### PART NUMBER

5N02-A300H